

Jellagen® 3D Scaffolds

Next Generation Jellyfish Collagen 3D Scaffolds for *in vitro* cell culture and tissue engineering.

PRODUCT DESCRIPTION

3D scaffolds of jellyfish collagen, suitable for cell culture research and tissue engineering purposes.

Product Numbers

- JSM96F, JSM96H, JSM96Q
- JSM48
- JSM24F, JSM24H
- JSM06

*bespoke and non cross-linked formats available on request and subject to volume

FEATURES AND BENEFITS

FEATURES	BENEFITS
Innovative	Offers a viable alternative to mammalian and synthetic scaffolds.
Non-mammalian	Highly purified jellyfish collagen alternative providing consistent, repeatable results.
Batch to batch consistency	Offers improved research productivity allowing security of product consistency and reproducible results.
Cross-linked	EDC cross-linked scaffolds for enhanced mechanical and thermal stability, resorbable <i>in vivo</i> ¹ .
Evolutionary ancient collagen demonstrating sequence homology to collagen I and II	Universal applications for multiple cell types and regenerative medicine.
Manufactured according to ISO13485	Follows a quality controlled manufacturing process producing a consistent scaffold.
Uniform pore size	Promotes cell seeding, invasion, proliferation and differentiation. Allows for growth factor permeation and gas exchange ensuring long-term cell survival.
Natural Scaffold	Jellagen® 3D Scaffolds exhibit similar physiological components and properties to the ECM of the <i>in vivo</i> micro-environment.

1. Jonathan P. Widdowson, Alex J. Picton, Valerie Vince, Chris J. Wright, Andrew Mearns-Spragg. "In vivo comparison of jellyfish and bovine collagen sponges as prototype medical devices". J Biomed Mater Res B Appl Biomater. 2018 May;106(4):1524-153

The grade of Jellagen® jellyfish collagen used to manufacture these scaffolds has been tested to verify its applicability for routine cell culture research using human primary and iPSC-derived cell lines. Jellagen® Jellyfish collagen has been shown to promote cellular attachment, proliferation and differentiation to develop functional matrices.

Cell lines that have been cultured successfully on Jellagen® jellyfish include, but are not limited to: Mesenchymal Stem Cells (MSC's), fibroblasts, hepatocytes, endothelial cells, keratinocytes, chondrogenic progenitor cells, Urine Derived Stem Cells (UDC's), cardiomyocytes, ovarian cancer cells, iPSC-derived microglia and HEK293T.

PRODUCT INFORMATION	JELLAGEN® 3D SCAFFOLDS
Format	6, 24, 48 and 96-well scaffolds, cast in plates
Collagen	Jellyfish collagen
Serum level	Serum free
Storage	Store at room temperature
Shelf life	Under evaluation
Plate polymer	Non-tissue culture treated, polystyrene and non-pyrogenic
Colour	White to off-white
Bioburden	Negative
Shipping conditions	Room temperature
pH	Approximately 7.0 – 7.4 when suspended in PBS or tissue culture media

Useful References

- Jonathan P. Widdowson, Alex J. Picton, Valerie Vince, Chris J. Wright, Andrew Mearns-Spragg. "In vivo comparison of jellyfish and bovine collagen sponges as prototype medical devices". *J Biomed Mater Res B Appl Biomater*. 2018 May;106(4):1524-153
- Sourour Addad, J.Exposito, C.Faye, S.Ricard-Blum, and C. Lethias. "Isolation, Characterization and Biological Evaluation of Jellyfish Collagen for Use in Biomedical Applications". *Marine Drugs*. 2011; 9(6): 967–983
- Xiaochen Cheng, Ziyu Shao, Chengbo Li, Lejun Yu, Mazhar Ali Raja, and Chenguang Liu "Isolation, Characterization and Evaluation of Collagen from Jellyfish *Rhopilema esculentum* Kishinouye for Use in Hemostatic Applications. *PLoS One*. 2017: 12 (1)
- Seiya Miura and Shigeru Kimura. "Jellyfish Mesogloea Collagen – characterisation of molecules AS $\alpha 1\alpha 2\alpha 3$ heterotrimers". *The Journal of Biological Chemistry*. 1985. Vol. 260, No. 28, Issue of December 5, pp. 15352-15356.
- Eun Song, So Yeon Kim, Taehoon Chun, Hyun-Jung Byun, Young Moo Lee. "Collagen scaffolds derived from a marine source and their biocompatibility". *Biomaterials* 27. 2006. 2951–2961
- Judith Sewing1, Matthias Klinger and Holger Notbohm. "Jellyfish collagen matrices conserve the chondrogenic phenotype in two- and three- dimensional collagen matrices.". *Journal of Tissue Engineering and Regenerative Medicine*. 2015 Research Article.
- Birgit Hoyer, Anne Bernhardt, Anja Lode, Sascha Heinemann, Judith Sewing, Matthias Klinger, Holger Notbohm, Michael Gelinsky." Jellyfish collagen scaffolds for cartilage tissue engineering." *Acta Biomaterialia* 10.2014. 883–892
- Marion Pugliano, Xavier Vanbellinthen, Pascale Schwinté, Nadia Benkirane-Jesseland Laetitia Keller. "Combined Jellyfish Collagen Type II, Human Stem Cells and Tgf- $\beta 3$ as a Therapeutic Implant for Cartilage Repair." *Journal of Stem Cell Research & Therapy*. 2017, 7:4
- Ayako Miki, Satomi Inaba, Takayuki Baba, Koji Kihira, Harumi Fukada and Masayuki Oda. "Structural and physical properties of collagen extracted from moon jellyfish under neutral pH conditions". *Bioscience, Biotechnology, and Biochemistry*, 2015 Vol. 79, No. 10, 1603–1607

DISCLAIMER

This product is for R&D use only and is not intended for human or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

FRM-88REV00

